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HSV 1/2 RQ Kit

Manual

For Real-Time PCR Detection and
Quantitation of Herpes Simplex
Type 1 and 2

For use with Rotor-Gene and StepOne
Research use only

NG-WI-ASL-06-200

Version 2.0

March 2017


Table of Contents:

1. Introduction	2
2. Kit Contents	3
3. Storage and Stability	3
4. General Precautions.....	3
5. Additionally Required Materials.....	4
6. Specimen, storage and transport.....	4
7. Interfering Substances	4
8. DNA isolation.....	5
9. Internal control (IC).....	5
10. Quantitation.....	5
11. PCR Protocol.....	6
12. Programming of the Rotor-Gene	6
13. Programming of StepOne(/Plus)	6
14. Programming other machines	7
15. Data Analysis: Rotor-Gene.....	9
16. Data Analysis: StepOne	10
17. Linear Range.....	11
18. Sensitivity.....	12

HSV 1/2 RQ kit is intended for the quantitative detection of HSV Type 1 and 2 DNA with Rotor-Gene or StepOne instrument. This kit is for research use only.

1. Introduction

Human Herpes Virus type 1 and 2 (HHV-1, HHV-2) known as Herpes simplex virus 1 and 2 (HSV-1, HSV-2) are large enveloped DNA viruses with a genome of about 150 kbp double stranded DNA. HSV is a ubiquitous agent and can cause a variety of infections which in adults are usually benign. Following active infection, it establishes a latent infection in the sensory and autonomic ganglia of the nervous system. It can be reactivated with several potential triggers like stress. If reactivated, retrograde transmission of the virus from the peripheral site to CNS along nerve axon can result in CNS infection which makes HSV the most common cause of viral encephalitis. In such cases, prompt laboratory diagnosis of HSV infection is essential for patient management and possible initiation of antiviral therapy and reducing the mortality rate. Quantitative PCR assays can estimate the extent of viral replication in CNS and would be helpful as a prognostic marker and in monitoring antiviral effects.

HSV 1/2 RQ kit provides a ready-to-use Real-Time PCR system for detection and quantitation of HSV-1 and HSV-2 with Rotor-Gene or StepOne machines. Currently Real-Time PCR provides the highest sensitivity and widest dynamic range among other methods for detection of HSV and viral load determination. In this method application of fluorescent probes allows detection of amplified product. Analysis of fluorescent kinetics also leads to quantification of the target sequence in the reaction without requiring post-amplification analysis, reducing the possibility of contamination with the PCR product.

The kit also incorporates an *Internal Control* (IC) to identify possible PCR inhibition.

2. Kit Contents

The kit contains a manual, a CD with Rotor-Gene and StepOne templates and following reagents:

Label	Content	Quantity
HSV Mix	PCR Master mix*	360 µl
HSV S1	Standard 1: 10,000 copy/µl	150 µl
HSV S2	Standard 2: 1,000 copy/µl	150 µl
HSV S3	Standard 3: 100 copy/µl	150 µl
HSV S4	Standard 4: 10 copy/µl	150 µl
Water	PCR Grade Water	200 µl

* 1, 2 and 4 tubes for 24, 48 and 96 reaction kits.

3. Storage and Stability

The kit components should be shipped and stored at -20°C and are stable until the expiry date mentioned. Avoid repeated freeze-thaws more than few times to prevent reduced sensitivity.

4. General Precautions

In order to prevent false results, always pay attention to the following points:

- Treat all samples as potentially infectious.
- Within the pre-PCR work area assign three separate spaces for: a) Sample storage and extraction, b) Reagent preparation where the HSV Mix is aliquoted into tubes and c) Reaction preparation area for addition of extracted DNA to the tubes.

- Always wipe the working surfaces with 70% Ethanol before and after work.
- Thaw kit components on **crushed ice** completely, mix by flickering followed by a quick spin and **store on crushed ice while working**.
- Do not place PCR tubes on crushed ice. Use cooling blocks instead.

5. **Additionally Required Materials**

To use this kit, you need the following items:

- Real-Time PCR machine and accessory computer
- Table top microtube centrifuge
- Vortex Mixer
- Dry Block Heater
- Adjustable pipettors and nuclease free filtered tips
- DNA extraction kit
- Nuclease free 1.7ml microtubes and PCR microtubes
- Disposable powder-free gloves
- Cold block

6. **Specimen, storage and transport**

CSF or swab of the genital lesions should be collected in sterile condition in proper and sterile tubes. Samples should be shipped at +4°C. Upon receipt sample can be stored at +4°C for few days or stored at -20°C for up to few weeks.

7. **Interfering Substances**

Heparin (more than 10 IU/ml) affects the PCR. Blood collected in heparin containing tubes should not be used. Samples of heparinized patients must not be used as well.

Elevated levels of bilirubin (≤ 4.5 mg/dl) and lipids (≤ 1000 mg/dl) and hemolytic samples do not influence the extraction and PCR.

8. Internal control (IC)

To examine the possible PCR inhibition and to prevent false negative results, an *Internal Control* is included in the HSV Mix. IC should generate a CT of 26-34 in Yellow Channel.

9. DNA isolation

DNA isolation can be performed using different kits from various manufacturers. We recommend the following:

- High Pure Viral Nucleic Acid Kit (Cat. no. 11858874001, Roche Applied Science, Mannheim, Germany)
- QIAamp DNA Blood Mini Kit (Cat. no. 51104, Qiagen GmbH, Hilden, Germany)
- QIAampUltraSens ® Virus Kit (Cat. no. 53704, Qiagen GmbH, Hilden, Germany)
- QIAampMiniElute Virus Spin Kit (Cat. no. 57704, Qiagen GmbH, Hilden, Germany).

10. PCR Protocol

Thaw the reagents on ice completely followed by a brief mixing and a quick spin. Place required number of tubes on cold block. Consider one tube for each sample plus one for each standard and one for the negative control.

Pipette 15ul of HSV Mix directly to each tube followed by adding 10ul of standard or isolated DNA.

Cap the tubes and visually inspect to make sure all are capped securely. Place tubes in the machine.

Note: Working with StepOne instrument, spin tubes briefly before loading on the block.

Note: If using Rotor-Gene attach the locking ring.

11. Programming Rotor-Gene

Before you start the machine, make sure you have attached the locking ring on the rotor!

Open the CD provided in the kit and double click on HSV 0.2 or HSV 0.1 strip according to the used tubes. Program starts. Click on Start button (Green button on the top menu). On the pop up window click start again and save the run file.

12. Programming StepOne

Open the StepOne software (V 2.*). On the Set Up menu click on Template and select the file on CD provided with the kit. Click on Plate Setup. One negative control, 4 standards and 10 samples are defined. You may change plate set up using right click options (copy, paste, clear). You may also add/remove samples or change sample names on “Define Targets and Samples” menu. When finished, click on “Start Run” and save the experiment on desired location. Instrument will start shortly.

13. Programming other machines

If you apply this kit to other Real-Time PCR machines, program it according to the following table:

Step	Temperature and time	Cycles
1	95°C x 10 min	1 cycle
2	95°C x 15 sec	45 cycles
	60°C x 60 sec	

Fluorescence should be collected at 60°C for FAM and VIC dyes. The HSV Mix contains ROX with the final concentration of 300nM in reaction.

14. Data Analysis: Rotor-Gene

Before analyzing results, make sure in the sample menu all the standards have been defined as "standard" and the relative concentrations have been entered. Patient samples should be defined as "unknown" and Negative control or no template control as "Negative Control" or "NTC" respectively.

Analyze data according to manufacturer recommendations. Perform quantitative analysis for **HSV (Green channel)** and qualitative analysis for **Internal Control (Yellow channel)**. Briefly, click on analysis menu and then under Quantitation tab double click on cycling A. Green.

In the pop up for Automatic Threshold increase the minimum or lower bound until it surpasses the negative control or the NTC fluorescence, and then click on OK.

Repeat the above for Cycling A. Yellow but, cancel the Automatic Threshold and manually put threshold on 0.1.

Figures 1 and 2 represent typical graphs for Rotor-Gene machine.


Figure 1: Typical HSV Graph in Green Channel for Rotor-Gene


Figure 2: Typical IC Graph in Yellow Channel for Rotor-Gene

Note that a sample is considered Positive only if it has a sigmoid graph and log phase, and only then CT is reliable and can be used.

In the absence of sigmoid graph and log phase, sample is considered Negative and CT if present is not reliable.

Consider following points when analyzing:

- A sample is **Positive** if it is positive in Green channel with a sigmoid graph and a CT of less than 40. The viral load or quantitation results in the Cycling A. Green are valid.
- A sample is **Negative** if it is negative in Green channel while it is positive in Yellow channel with a sigmoid graph and CT of 26-34.
- Results are **Inconclusive** and the test should be repeated if a sample is negative in both of Green and Yellow channels.

15. Data Analysis: StepOne

Analyze data according to manufacturer recommendations. Briefly, click on Analyze and set the threshold for **HSV/FAM** at 0.1 and at 0.05 for **IC/VIC**.

Figures 3 and 4 represent typical graphs for StepOne machine.


Figure 3: Typical HSV Graph in FAM Channel for StepOne


Figure 4: Typical IC Graph in VIC Channel for

Note that a sample is considered Positive only if it has a sigmoid graph and log phase, and only then CT is reliable and can be used.

In the absence of sigmoid graph and log phase, sample is considered Negative and CT if present is not reliable.

Consider following points when analyzing:

- A sample is **Positive** if it is positive in FAM/HSV channel with a sigmoid graph and a CT of less than 40. The viral load or quantitation results are valid.
- A sample is **Negative** if it is negative in FAM/HSV channel while it is positive in VIC/IC channel with a sigmoid graph and CT of 26-34.
- Results are **Inconclusive** and the test should be repeated if a sample is negative in both of FAM/HSV and VIC/IC channels.

16. Quantitation

The kit provides 4 quantitation standards with defined titers to generate a standard curve for quantification of samples viral load. Working with Rotor-Gene machine, the standard curve from a previous run can also be imported for quantification of samples to the recent run. To do so, at least one standard must be used in the current run. Apparently using all five standards in each run will lead to more accurate results.

Quantitation standards are defined as copy/ul. To convert the result to copy/ml following equation should be used:

$$\text{Result(copy/ml)} = \frac{\text{Result(copy/}\mu\text{l)} \times \text{elution volume}(\mu\text{l})}{\text{sample volume(ml)}}$$

“Sample volume” is the sample volume used for DNA isolation and “Elution volume” is the volume of buffer or water used to elute or dissolve isolated DNA.

17. Linear Range

The linear range of the kit was assessed with dilution series of the cloned target and showed to be linear in the range of 1,000,000 copy/ μ l to 2 copies/ μ l.

18. Sensitivity

The analytical detection limit of the kit was assessed with dilution series of the cloned target and showed a limit of detection equal to 0.5 copy/ μ l.