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# **VZV RQ Kit Manual**

**For Real-Time PCR Detection and  
Quantitation of Varicella Zoster virus  
(VZV) DNA**

For use with Rotor-Gene or StepOne  
Research use only

NG-WI-ASL-41  
Version 1.0  
Spring 2021

## Table of Contents:

|                                          |    |
|------------------------------------------|----|
| 1. Introduction .....                    | 2  |
| 2. Kit Contents .....                    | 2  |
| 3. Storage and Stability .....           | 3  |
| 4. General Precautions.....              | 3  |
| 5. Additionally Required Materials ..... | 4  |
| 6. Specimen, storage and transport.....  | 4  |
| 7. Interfering substances .....          | 4  |
| 8. Internal control (IC).....            | 5  |
| 9. DNA isolation.....                    | 5  |
| 10. PCR Protocol.....                    | 6  |
| 11. Programming of the Rotor-Gene .....  | 6  |
| 12. Programming of StepOne(/Plus) .....  | 7  |
| 13. Programming other machines .....     | 7  |
| 14. Data Analysis: Rotor-Gene.....       | 8  |
| 15. Data Analysis: StepOne .....         | 10 |
| 16. Quantitation.....                    | 12 |
| 17. Linear Range.....                    | 13 |
| 18. Sensitivity.....                     | 13 |

**VZV RQ kit** is intended for use with Real-Time PCR machine and for the quantitative detection of Varicella Zoster virus (VZV) DNA. This kit is for research use only.

## 1. Introduction

Varicella Zoster virus (VZV) is a double-stranded 125 kb DNA virus and a member of Herpesviridae family. Varicella (chickenpox) is the result of childhood infection. Following the initial infection virus remains latent in nervous ganglia. Reactivation of the virus in adults results in Zoster (shingles). Infection during pregnancy may also lead to congenital defects.

VZV RQ kit provides a ready-to-use Real-Time PCR system for detection and quantitation of Varicella Zoster virus with Rotor-Gene or StepOne machines. This method provides the highest sensitivity and widest dynamic range among other available methods. In this method application of fluorescent dye labeled probes allows detection of amplified product. Analysis of fluorescent kinetics also leads to quantification of the target sequence in the reaction without requiring post-amplification analysis, reducing the possibility of contamination with the PCR product.

This kit also incorporates an *Internal Control* (IC) to identify possible extraction failure or PCR inhibition.

## 2. Kit Contents

The kit contains a manual, a CD with Rotor-Gene and StepOne templates and following reagents:

| Label   | Content                    | Quantity |
|---------|----------------------------|----------|
| VZV Mix | PCR Master mix*            | 360 µl   |
| VZV S1  | Standard 1: 10,000 copy/µl | 150 µl   |
| VZV S2  | Standard 2: 1,000 copy/µl  | 150 µl   |
| VZV S3  | Standard 3: 100 copy/µl    | 150 µl   |
| VZV S4  | Standard 4: 10 copy/µl     | 150 µl   |
| VZV IC  | Internal control           | 250 µl   |
| Water   | PCR Grade Water            | 200 µl   |

\* 1, 2 and 4 tubes for 24, 48 and 96 reaction kits, respectively.

### 3. Storage and Stability

The kit components should be shipped and stored at -20°C and are stable until the expiry date mentioned. Avoid repeated freeze-thaws more than few times to prevent reduced sensitivity.

### 4. General Precautions

In order to prevent false results, always pay attention to the following points:

- Treat all samples as potentially infectious.
- Within the pre-PCR work area assign three separate spaces for: a) Sample storage and extraction, b) Reagent preparation where the master-mix is aliquoted into tubes and c) Reaction preparation area for addition of extracted DNA to the tubes.
- Always wipe the working surfaces with 70% Ethanol before and after work.

- Thaw kit components **on ice** completely, mix by flickering followed by a quick spin and **store on crushed ice** after.
- Do not place PCR tubes on crushed ice. Use cooling blocks instead.

## 5. Additionally Required Materials

To use this kit, you need the following items:

- Real-Time PCR machine and accessory computer
- Table top microtube centrifuge
- Vortex Mixer
- Dry Block Heater
- Adjustable pipettors and nuclease free filtered tips
- DNA extraction kit
- Nuclease free 1.7ml microtubes and PCR microtubes
- Disposable powder-free gloves

## 6. Specimen, Storage and Transport

Peripheral blood plasma, CSF or tissue should be collected in sterile condition in proper and sterile containers. Samples can be stored at +4°C for few days or stored at -20°C for up to few weeks. The minimum sample volume recommended for testing is 200µl of plasma, which requires 0.5ml of whole blood.

## 7. Interfering Substances

Heparin (more than 10 IU/ml) affects the PCR. Blood collected in heparin containing tubes should not be used. Samples of heparinized patients must not be used as well. Elevated levels of

bilirubin ( $\leq 4.5$  mg/dl) and lipids ( $\leq 1000$  mg/dl) and hemolytic samples do not influence the extraction and PCR.

## 8. Internal Control (IC)

In order to evaluate the possibility of DNA extraction failure and PCR inhibition and prevent false negative results, VZV RQ kit contains internal control (IC). IC can be used during extraction process or simply added to VZV Mix. To monitor both DNA extraction and PCR reaction, internal control should be added to the mixture of lysis buffer and patient sample. Required volume of internal control is 10% of elution buffer. For instance if extracted DNA is eluted with 100ul, then 10ul of internal control should be added to the mixture of lysis buffer and patient sample.

**Please note that internal control should not be added directly to the patient sample (i.e. before addition of lysis buffer) as it loses its efficiency.**

If IC is added to VZV Mix, only PCR inhibition can be monitored. For this purpose, 1ul of IC reagent should be added to each reaction. For example, for 10 PCR reactions, 15ul of internal control should be added to 150ul of VZV Mix before it is added to the tubes. In a successful DNA extraction and PCR test, Internal control should generate a CT of 28-32 in Yellow Channel on Rotor-Gene and a CT of 28-34 in VIC channel on StepOne.

## 9. DNA Isolation

DNA isolation can be performed using different kits from various manufacturers. We recommend the following:

- High Pure Viral Nucleic Acid Kit (Cat. no. 11858874001, Roche Applied Science, Mannheim, Germany)
- QIAamp DNA Blood Mini Kit (Cat. no. 51104, Qiagen GmbH, Hilden, Germany)
- QIAampUltraSens ® Virus Kit (Cat. no. 53704, Qiagen GmbH, Hilden, Germany)
- QIAampMiniElute Virus Spin Kit (Cat. no. 57704, Qiagen GmbH, Hilden, Germany).

## 10. PCR Protocol

Thaw the reagents on ice completely followed by a brief mixing and a quick spin. Place required number of tubes on cold block. Consider one tube for each sample plus one for each standard and one for the negative control.

**If IC is introduced during extraction process, pipette 15µl of VZV Mix to each PCR tube.**

**If IC is added to the VZV Mix, add 15ul of the prepared mix (as described in section 8) to each PCR tube.**

**Then add 10ul of extracted DNA, standard or water to each tube.**

Cap the tubes and visually inspect to make sure all are capped securely. Place tubes in the machine.

*Note: Working with StepOne instrument, spin tubes briefly before loading on the block.*

*Note: If using Rotor-Gene attach the locking ring.*

## 11. Programming Rotor-Gene

*Before you start the machine, make sure you have attached the locking ring on the rotor!*

Open the CD provided in the kit and double click on VZV Template 0.1 or VZV 0.2 according to the tubes used. Program starts. Click on Start button (Green button on the top menu). On the pop up window click start again and save program on desired location.

## 12. Programming StepOne

Open the StepOne software (V 2.\*). On the Set Up menu click on Template and select the file on CD provided with the kit. Click on Plate Setup. One negative control, 4 standards and 10 samples are defined. You may change plate set up using right click options (copy, paste, clear). You may also add/remove samples or change sample names on “Define Targets and Samples” menu. When finished, click on “Start Run” and save the experiment on desired location. Instrument will start shortly.

## 13. Programming Other Machines

If you apply this kit to other Real-Time PCR machines, program it according to the following table:

| Step | Temperature and time | Cycles |
|------|----------------------|--------|
| 1    | <b>95°C x 10 min</b> | 1      |
| 2    | <b>95°C x 15 sec</b> | 45     |
|      | <b>60°C x 60 sec</b> |        |

Fluorescence should be collected at 60°C for FAM and VIC dyes. The VZV Mix contains ROX. Final concentration of ROX in reaction is 300nM.

#### 14. Data Analysis: Rotor-Gene

Before analyzing results, make sure in the sample menu all the standards have been defined as "standard" and the relative concentrations have been entered. Patient samples should be defined as "unknown" and Negative control or no template control as "Negative Control" or "NTC" respectively.

Analyze data according to manufacturer recommendations. Perform quantitative analysis for **VZV (Green channel)** and qualitative analysis for **Internal Control (Yellow channel)**. Briefly, click on analysis menu and then under Quantitation tab double click on cycling A. Green.

In the pop up for Automatic Threshold increase the minimum or lower bound until it surpasses the negative control or the NTC fluorescence, and then click on OK.

Repeat the above for Cycling A. Yellow but cancel the Automatic Threshold and manually put threshold on 0.1. Figures 1 and 2 represent typical graphs for Rotor-Gene machine.

**Note that a sample is considered Positive only if it has a sigmoid graph and log phase, and only then CT is reliable and can be used. In the absence of sigmoid graph and log phase, sample is considered Negative and CT if present is not reliable.**

Consider following points when analyzing:

- A sample is **Positive** if it is positive in Green channel with a sigmoid graph and a CT of less than 40. The viral load or quantitation results in the Cycling A. Green are valid.
- A sample is **Negative** if it is negative in Green channel while it is positive in Yellow channel with a sigmoid graph and CT of 28-32.
- Results are **Inconclusive** and the test should be repeated if a sample is negative in both of Green and Yellow channels.


Fig 1. Typical VZV graph in Green channel for Rotor-Gene


Fig 2. Typical IC graph in Yellow channel for Rotor-Gene


## 15. Data Analysis: StepOne

Analyze data according to manufacturer recommendations. Briefly, click on Analyze and set the threshold for VZV/FAM and IC/VIC at 0.1. Figures 3 and 4 represent typical graphs for StepOne machine.


**Note that a sample is considered Positive only if it has a sigmoid graph and log phase, and only then CT is reliable and can be used. In the absence of sigmoid graph and log phase, sample is considered Negative and CT if present is not reliable.**

Consider following points when analyzing:

- A sample is **Positive** if it is positive in FAM/VZV channel with a sigmoid graph and a CT of less than 40. The viral load or quantitation results are valid.
- A sample is **Negative** if it is negative in FAM/VZV channel while it is positive in VIC/IC channel with a sigmoid graph and CT of 28-34.
- Results are **Inconclusive** and the test should be repeated
- if a sample is negative in both of FAM/VZV and VIC/IC channels.

**Fig 3.** Typical VZV graph in FAM channel for StepOne

**Fig 4.** Typical IC graph in VIC channel for StepOne

## 16. Quantitation

The kit provides 4 quantitation standards with defined titers to generate a standard curve for quantification of samples viral load. Working with Rotor-Gene machine, the standard curve from a previous run can also be imported for quantification of samples to the recent run. To do so, at least one standard must be used in the current run. Apparently using all five standards in each run will lead to more accurate results.

Quantitation standards are defined as copy/ $\mu$ l. To convert the result to copy/ml following equation should be used:

$$\text{Result (copy/ml)} = \frac{\text{Result(copy/}\mu\text{l)} \times \text{elution volume (}\mu\text{l)}}{\text{sample volume (ml)}}$$

“Sample volume” is the plasma volume used for DNA isolation and  
“Elution volume” is the volume of buffer or water used to elute or  
dissolve isolated DNA.

## **17. Linear Range**

The linear range of the kit was assessed with dilution series of the cloned target and showed to be linear in the range of 100,000,000 copies/ $\mu\text{l}$  to 10 copy/ $\mu\text{l}$ .

## **18. Sensitivity**

The analytical detection limit of the kit was assessed with dilution series of the cloned target and showed a limit of detection equal to 4 copy/ $\mu\text{l}$ .