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PAI-1 RQ Kit Manual

**For Real-Time PCR Detection of PAI-1 4G/5G
polymorphism**

For use with StepOne and Rotor-Gene
Research use only

NG-WI-ASL-22-101

Version 1.2

March 2019


Table of Contents:

1. Introduction	2
2. Kit Contents.....	3
3. Storage and Stability	3
4. General Precautions.....	3
5. Additionally Required Materials	4
6. Specimen, Storage and Transport	4
7. Interfering Substances.....	4
8. Internal Control	4
9. DNA Isolation	5
10. PCR Protocol.....	5
11. Programming Rotor-Gene.....	5
12. Programming StepOne	5
13. Programming Other Machines.....	6
14. Data Analysis: Rotor-Gene.....	6
15. Data Analysis: StepOne	8

PAI-1 RQ kit is for the detection of 4G/5G polymorphism in promoter region of PAI-1 gene in human genomic DNA. This kit is intended for research use only and is designed for StepOne or Rotor-Gene machine.

1. Introduction

PAI-1 (Plasminogen Activator Inhibitor 1) is a 52 kDa glycoprotein. It is the main inhibitor of activators of plasminogen, and hence fibrinolysis and blood clot degradation. PAI-1 gene is located on chromosome 7 and consists of 9 exons. It includes many variations. The most significant variation, however, lies in the promoter region at -675, which is a single guanine base insertion/deletion known as 4G/5G polymorphism. While 5G allele binds both transcription activator and repressor, 4G allele only binds activator and leads to higher expression of PAI-1. 4G allele is considered a predisposing factor and has been related to various clinically important conditions including deep vein thrombosis, atherosclerosis and myocardial infarction, cancer, sepsis, miscarriage and asthma.

PAI-1 RQ kit provides a ready-to-use Real-Time PCR system for detection of PAI-1 4G/5G polymorphism with StepOne or Rotor-Gene machine.

This method applies fluorescent probes allowing detection of amplified product. Analysis of fluorescent kinetics also leads to assessment of the target sequence in the reaction without requiring post-amplification analysis, reducing the possibility of contamination with the PCR product.

2. Kit Contents

The kit contains a manual, a CD with StepOne templates and following reagents:

Label	Content	Quantity
PAI-1 RQ Mix	Ready to use PCR Mix*	480 µl
4G/4G Ctrl	4G/4G Homozygous Control	50 µl
4G/5G Ctrl	4G/5G Heterozygous Control	50 µl
5G/5G Ctrl	5G/5G Homozygous Control	50 µl
Water	PCR Grade Water	200 µl

* 1, 2 or 4 tubes for 24, 48 or 96 reaction kits, respectively.

3. Storage and Stability

The kit components should be shipped and stored at -20°C and are stable until the expiry date mentioned. Avoid repeated freeze-thaws more than few times to prevent reduced sensitivity.

4. General Precautions

In order to prevent false results, always pay attention to the following points:

- Treat all samples as potentially infectious.
- Within the pre-PCR work area assign three separate spaces for: a) Sample storage and extraction, b) Reagent preparation where the master-mix is aliquoted into tubes and c) Reaction preparation area for addition of extracted DNA to the tubes.
- Always wipe the working surfaces with 70% Ethanol before and after work.
- Thaw kit components on crushed ice completely, mix by flickering followed by a quick spin and store on crushed ice after.
- Do not place 0.2ml PCR tubes on crushed ice. Use cold blocks instead.

5. Additionally Required Materials

To use this kit, you need the following items:

- Real-Time PCR machine and accessory computer
- Table top microtube centrifuge
- Vortex Mixer
- Dry Block Heater
- Adjustable pipettors and nuclease free filtered tips
- DNA extraction kit
- Nuclease free 1.7ml microtubes and PCR microtubes
- Disposable powder-free gloves
- Cold block

6. Specimen, Storage and Transport

Required sample for this test is peripheral blood. Sample should be collected in sterile condition in proper tubes containing EDTA or citrate as anticoagulant. Whole blood should be shipped and stored at +4°C (stable for 3 days). For long term storage, sample should be aliquoted and stored at -20°C which remains stable for few months. DNA can be extracted from whole blood or buffy coat.

7. Interfering Substances

Heparin (more than 10 IU/ml) affects the PCR. Blood collected in heparin containing tubes should not be used. Samples of heparinized patients must not be used as well. Elevated levels of bilirubin (≤ 4.5 mg/dl) and lipids (≤ 1000 mg/dl) and hemolytic samples do not influence the extraction and PCR.

8. Internal Control (IC)

Either of wild type or mutant alleles may function as internal control for this assay as any patients carry 4G or 5G or both alleles. If a patient is negative for both, then the test is not valid and should be

repeated.

9. DNA Isolation

DNA isolation can be performed using different kits from various manufacturers. We recommend the following:

- High Pure PCR Template Preparation Kit (Cat# 11796828001, Roche Applied Science, Mannheim, Germany)
- QIAamp DNA Blood Mini Kit (Cat. no. 51104, Qiagen GmbH, Hilden, Germany)

10. PCR Protocol

Thaw the reagents on ice completely followed by a brief mixing and a quick spin. Place required number of tubes on cold block. Consider one tube for each sample plus four for controls and NTC.

Please note that the software can call the genotype of unknown samples only in the presence of all the controls provided with the kit.

Pipette 20µl of PAI-1 RQ Mix directly to each tube followed by adding 5µl of sample or control.

Cap the tubes and visually inspect to make sure all are capped securely. Place tubes in the machine.

Note: Working with StepOne instrument, spin tubes briefly before loading on the block.

11. Programming Rotor-Gene

- Before you start the machine, make sure you have attached the locking ring on the rotor!

In the provided CD double click on PAI-1 RQ 0.1 or PAI-1 RQ 0.2 (according to the used microtubes) to open the program. Click on Start button (Green button on the top menu). On the pop up window click start again and save the run file.

Edit sample names.

12. Programming StepOne

Open the StepOne software (V 2.*). On the Set Up menu click on Template and select the file on CD provided with the kit. Click on Plate Setup. One negative control, 3 positive controls and 10 samples are defined. You may change plate set up by using right click options (copy, paste, clear). You may also add /remove samples or change sample name on “Define Targets and Samples” menu. When finished, click on “Start Run” and save the experiment on desired location. Instrument will start shortly.

13. Programming Other Machines

If you apply this kit to other Real-Time PCR machines, program it according to the following table:

Step	Temperature and time	Cycles
1	95°C x 10 min	1
2	95°C x 15 sec	50
	65°C x 40 sec	

Fluorescence should be collected at 65°C for FAM and VIC dyes. PAI-1 RQ Mix contains 300nM ROX in final reaction.

14. Data Analysis: Rotor-Gene

Analyze data according to manufacturer recommendations. Briefly, click on analysis menu select “Other” and then “Scatter Graph Analysis”. Using Ctrl button, mark both Green and Yellow channels and then click on “Show”.

The “Scatter Graph Window” is now shown in the lower left window. Each dot represents one sample. The vertical axis shows the

Green fluorescence and horizontal axis is for the Yellow fluorescence. Note that 4G allele is detected in Green channel and 5G allele in Yellow channel.

Therefore, samples are located in four regions of the graph. 4G/4G samples with high Green and low Yellow fluorescence gather in upper left; 5G/5g samples with low Green fluorescence and high Yellow gather in lower right; 4G/5G samples with average Green and Yellow fluorescence gather in upper right; and finally NTC or control with no DNA gather in lower left with low Green and low Yellow fluorescence.


Figure 1: Defining genotypes of different controls of the test on Rotor-Gene instrument.

In order to determine genotypes of the samples, above regions should be defined on the graph. To do so, turn off all the samples except for the negative and positive controls. Then, holding the right click, draw a rectangle around each control. Each box

represents one of the above mentioned regions. Label 4G/4G genotype region as “Mutant”, 4G/5G as “Heterozygous”, 5G/5G as “Wild Type” and Negative Control as “None”(Figure 1).

15. Data Analysis: StepOne

Analyze data according to manufacturer recommendations. Briefly, click on Analyze and select “Allelic Discrimination” tab. Genotype of the samples are automatically determined by the software based on comparison of their fluorescence with the controls.

Please note that the software can call the genotype of unknown samples only in the presence of all the controls provided with the kit!

Figure 3 represents the typical graph. Each dot represents one sample. All samples are color coded based on their genotype as blue for 4G/4G, Red for 5G/5G, Green for 5G/4G, and Black box for control with no DNA or NTC. Unknown samples will be marked with black "x".

Note that the vertical axis shows the FAM fluorescence and horizontal axis is for the VIC fluorescence. Accordingly, 4G allele is detected in FAM channel and 5G allele in VIC channel. Therefore, samples are located in four regions of the graph. 4G/4G samples with high FAM and low VIC fluorescence gather in upper left; 5G/5G samples with low FAM fluorescence and high VIC gather in lower right; 5G/4G samples with high FAM and high VIC fluorescence gather in upper right; and finally NTC or control with no DNA gather in lower left with low FAM and low VIC fluorescence (Figure 2).

Note! If location of controls are not similar to the typical graph represented here; or if controls are located too close to each other in a graph; and if samples lie in between controls and genotypes cannot be called, results are invalid and test should be repeated!


Figure 2: Scatter graph of positive and negative controls on StepOne.