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MTHFR-677 RG

Kit Manual

**For Real-Time PCR Detection of MTHFR
(C677T) Variant**

For use with Rotor-Gene 6000/Q
Research use only

NG-WI-ASL-20-101
Version 1.2
December 2019


Table of Contents:

1. Introduction	2
2. Kit Contents	3
3. Storage and Stability	3
4. Additionally Required Materials	3
5. General Precautions	4
6. Specimen, Storage and Transport	5
7. Interfering Substances	5
8. DNA Isolation	5
9. Internal Control (IC)	6
10. PCR Protocol	6
11. Programming of RotorGene6000/Q	7
12. Data Analysis	7

MTHFR-677 RG kit is intended for the detection of MTHFR (C677T) variant in genomic DNA extracted from human cells with RotorGene machine. This kit is intended for research use only

1. Introduction

MTHFR stands for Methylenetetrahydrofolate reductase which is a key regulatory enzyme in folate and homocysteine metabolism. Since, the C677T and A1298C variants show reduced activity (30-65% and 15-30% respectively), these polymorphisms are among the the most important genetic determinants of plasma total homocysteine. Due to their influence on homocysteine metabolism, they have been linked to several complex disorders including arterial thrombosis, atherosclerosis, coronary heart disease, neural tube defects, several types of cancer, adverse drug interactions and osteoporosis though, findings are still controversial. Depending on the ethnic background up to 50% of population may carry the variant allele of MTHFR.

MTHFR-677 RG kit provides a ready-to-use Real-Time PCR system for detection of C677T variant with RotorGene machine.

This method applies fluorescent probes allowing detection of amplified product. Analysis of fluorescent kinetics also leads to assessment of the target sequence in the reaction without requiring post-amplification analysis, reducing the possibility of contamination with the PCR product.

2. Kit Contents

The kit contains a manual, a CD with Rotor-Gene™ template and following reagents:

Label	Content	Quantity
MTHFR Mix	PCR Master mix*	480µl
MM Pos Ctrl	MM Positive Control	50µl
MW Pos Ctrl	MW Positive Control	50µl
WW Neg Ctrl	WW Negative Control	50µl
Water	PCR Grade Water	200µl

* 1, 2 or 4 tubes for 24, 48, or 96 reaction kits, respectively.

3. Storage and Stability

All the kit components should be shipped and stored at -20°C and are stable until the expiry date mentioned. Avoid repeated freeze-thaws more than 3 times to prevent reduced sensitivity.

4. Additionally Required Materials

To use this kit, you need the following items:

- Real-Time PCR machine and accessory computer
- Table top microtube centrifuge
- Vortex Mixer
- Dry Block Heater
- Adjustable pipetters and nuclease free filtered tips
- DNA extraction kit
- Nuclease free 1.7ml microtubes and PCR microtubes
- Disposable powder-free gloves
- Cold block

5. General Precautions

In order to prevent false results, always pay attention to the following points:

- **Treat all samples as potentially infectious.**
- Within the pre-PCR work area assign three separate spaces for a) Sample storage and extraction; b) Reagent preparation where the Master Mix is aliquoted into tubes; and c) Reaction preparation area for addition of extracted samples to the tubes.
- Always wipe the working surfaces with 70% Ethanol before and after work.

- Thaw kit components on **crushed ice** completely, mix by flickering followed by a quick spin and **store on crushed ice** after.
- Do not place PCR tubes on crushed ice. Use cold blocks instead.

6. Specimen, Storage and Transport

Peripheral blood should be collected in sterile condition in proper and sterile tubes. We recommend EDTA or citrate as anticoagulant for MTHFR C677T variant detection. Whole blood should be shipped and stored at +4°C (stable for 3 days). Whole blood or buffy coat can be stored at +4°C for a few days. Otherwise should be aliquoted and stored at -20°C which is stable for few months.

DNA can directly be extracted from whole blood or buffy coat.

7. Interfering Substances

Heparin (more than 10 IU/ml) affects the PCR. Blood collected in heparin containing tubes should not be used. Samples of heparinized patients must not be used as well. Elevated levels

of bilirubin (≤ 4.5 mg/dl) and lipids (≤ 1000 mg/dl) and hemolytic samples do not influence the extraction and PCR.

8. DNA Isolation

DNA isolation can be performed using different kits from various manufacturers. We recommend using either High Pure PCR Template Preparation Kit (Cat. no. 11796828001, Roche Applied Science, Mannheim, Germany) or QIAamp DNA Blood Mini Kit (Cat. no. 51104, Qiagen GmbH, Hilden, Germany)

9. Internal Control (IC)

This kit requires no internal control since all the patients carry either wild type or mutant copy of the *MTHFR* gene. Therefore, it must be positive for at least wild type or the mutant copy of the gene. If a patient is negative for both above, then the test should be repeated.

10. PCR Protocol

Thaw the reagents on ice completely followed by brief mixing and a quick spin. Place required number of tubes on cold block.

MTHFR-677 RG (V1.1)

Consider one tube for each sample plus one for each control and one for no-template-control (NTC).

Pipette 20ul of MTHFR Mix directly to each tube followed by adding 5ul of standard or isolated DNA.

Cap the tubes and visually inspect to make sure all are capped securely. Place tubes in the machine.

Note: If using Rotor-Gene attach the locking ring too.

11. Programming of RotorGene 6000/Q

- Before you start the machine make sure you have attached the ring on the rotor!

- Please note that at temperatures of 30°C and above, RotorGene cooling system fails and test results are not valid! The optimum temperature for proper function of this machine is 20-25°C!

Open the CD provided in the kit and double click on MTHFR 677 RG 0.2 or on MTHFR 677 Strip (according to the used tubes) to open the program. Click on Start button (Green button on the top menu). On the pop up window click start again and save program. On the pop up sample window, name the samples according to the list.

You may also set Rotor-Gene according to the following table:

Step	Temperature and time	Cycles
------	----------------------	--------

1	95°C x 3 min	1
2	95°C x 15 sec	40
	55°C x 40 sec*	
	72°C x 20 sec	
3	95°C x 15 sec	1
4	40°C x 10 min	1
5	45°C x 1 min	1
6	Melt 50°C-80°C *	1

Fluorescence should be collected at steps marked with star (*) for FAM.

The PCR Mix is ROX free.

12. Data Analysis

Analyze the data according to manufacturer recommendations. Briefly, on analysis menu click on Melt tab and then, double click on “Melt A Green”.

On the right side menu, uncheck the box for “Flip sign of dF/dT”. On the bottom of “Melt Curve Analysis” window click on “Auto-Scale” tab. Set threshold manually above the background, so only two picks remain above the threshold. **The left pick (at 60-61°C) is for wild type allele (W) and the right one (at 66-67°C) is for the mutant allele (M).** Analyze the results as follows:

- A sample with only a single pick at 60-61°C is Negative for C677T variant and is considered **Negative or wild type homozygous (WW)**.
- A sample with only a single pick at 66-67°C is **Positive and homozygous (MM) for C677T** variant (both alleles are C677T variant).
- A sample with two picks at 60-61°C and at 66-67°C is **Positive for C677T variant and heterozygous (WM)**. Such a patient carries both wild type and variant alleles for MTHFR.

Results are inconclusive and the test should be repeated if:

- No pick is detected for a patient at either 60-61 or 66-67°C.
- A pick is detected for NTC.

Genotypes can be called automatically if bins are defined for W and M picks. Please refer to the RotorGene 6000/Q manual for further information.

Figures 1, 2 and 3 represent sample analysis and the relevant peak for each genotype.


Figure 1: Wild type (WW) graph on RotorGene


Figure 2: Mutant Heterozygote (WM) graph on RotorGene


Figure 3: Mutant Homozygote (MM) graph on RotorGene