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Flu-RSV RQ Kit Manual

**For Real-Time RT-PCR Detection
of Influenza A, Influenza B and RSV**

For use with Rotor-Gene
Research use only

NG-WI-ASL-39-100
Version 1.0
Autumn 2020


Table of Contents:

1. Introduction	2
2. Kit Contents	2
3. Storage and Stability	3
4. General Precautions.....	3
5. Additionally Required Materials	4
6. Specimen, storage and transport.....	4
7. Internal Control (IC)	4
8. RNA isolation.....	4
9. RT-PCR Protocol	5
10. Programming of the Rotor-Gene	5
11. Programming other machines	6
12. Data Analysis: Rotor-Gene.....	6
13. Sensitivity.....	10

Flu-RSV RQ kit is intended for detection of RNA of Influenza A and B viruses as well as RSV. It is designed for use with Rotor-Gene. This kit is for research use only.

1. Introduction

Human influenza A and B are enveloped RNA viruses of the Orthomyxoviridae family. While both of these viruses cause seasonal flu epidemics, Influenza A virus is responsible for majority of them and may also cause pandemics.

Respiratory Syncytial Virus (RSV) is also an enveloped RNA virus and belongs to Paramyxoviridae family. RSV mostly causes lower respiratory tract infections with cold like symptoms in healthy individuals but it may cause also severe respiratory infections in some.

Flu-RSV RQ kit provides a ready-to-use One-Step Multiplex Real-Time RT-PCR system for detection of Influenza A virus, Influenza B virus and RSV with Rotor-Gene. In this method application of fluorescent probes allows detection of multiple targets within the same PCR reaction. Analysis of fluorescent kinetics also leads to detection of target sequences in the reaction without requiring post-amplification analysis therefore, reducing the possibility of contamination with the PCR product.

2. Kit Contents

The kit contains a manual a CD with templates and following reagents:

Label	Content	Quantity
Flu-RSV Mix	RT-PCR Master mix	360 μ l*
Pos Ctrl	Positive Control	100 μ l
Water	PCR Grade Water	200 μ l

* 1, 2 and 4 tubes for 24, 48 and 96 reaction kits respectively.

3. Storage and Stability

The kit components should be shipped and stored at -20°C and are stable until the expiry date mentioned. Avoid repeated freeze-thaws especially for Flu-RSV Mix more than few times to prevent reduced sensitivity.

4. General Precautions

In order to prevent false results, always pay attention to the following points:

- Treat all samples as potentially infectious.
- Within the pre-PCR work area assign three separate spaces for: a) Sample storage and extraction, b) Reagent preparation where the master-mix is aliquoted into tubes and c) Reaction preparation area for addition of extracted RNA to the tubes.
- Always wipe the working surfaces with 70% Ethanol before and after work.
- Thaw on ice kit components completely, mix by flickering followed by a quick spin and store on crushed ice after.
- Keep RT-PCR Mix tube at -20°C at all times. Take it out just before use and return it to freezer immediately after.
- Do not place 0.2ml PCR tubes on crushed ice. Use cold block instead.

5. Additionally Required Materials

To use this kit, you need the following items:

- Real-Time PCR machine and accessory computer
- Vortex Mixer
- Adjustable pipettors and nuclease free filtered tips
- RNA extraction kit and required equipments/items
- PCR microtubes
- Disposable powder-free gloves
- Cold block

This kit does not require cDNA synthesis reagents!

6. Specimen, Storage and Transport

We recommend upper respiratory tract samples mainly nasopharyngeal swabs. Samples can be stored at 2-8°C for few hours or at -20°C or lower for up to few days.

7. Internal Control (IC)

In order to evaluate the possibility of RNA extraction failure and PCR inhibition and prevent false negative results, Flu-RSV Mix contains primers and probe of a housekeeping gene (RNase P) as an internal control. In a successful RNA extraction and RT-PCR test, IC should generate a CT of 20-35 in Orange/Rox Channel.

8. RNA Isolation

RNA isolation can be performed using different kits from various manufacturers. We recommend the following:

- High Pure Viral Nucleic Acid Kit (Cat. no. 11858874001, Roche Applied Science, Mannheim, Germany)
- QIAamp Viral RNA Mini Kit (Cat. no. 52904, Qiagen GmbH, Hilden, Germany)

- QIAampUltraSens ® Virus Kit (Cat. no. 53704, Qiagen GmbH, Hilden, Germany)
- QIAampMiniElute Virus Spin Kit (Cat. no. 57704, Qiagen GmbH, Hilden, Germany).

9. RT-PCR Protocol

Thaw the reagents on ice completely followed by a brief mixing and a quick spin. Place required number of tubes on cold block. Consider one tube for each sample plus one for positive sample and one for the negative control/water.

Aliquot 15µl of Flu-RSV Mix directly to each PCR tube. Then add 5ul of extracted RNA, Positive Control or water.

Lower reaction volumes may also be used. Including:

Cap the tubes and visually inspect to make sure all are capped securely. Place tubes in the machine.

Note: If using Rotor-Gene attach the locking ring too.

10. Programming of the Rotor-Gene

Before you start the machine, make sure you have attached the locking ring on the rotor!

Open the CD provided in the kit and double click on "Flu-RSV 0.2" or on "Flu-RSV 0.1" depending on the tubes used.

Program starts. Click on Start button (Green button on the top menu). On the pop up window click start again and save program on desired location. You may also set the machine according to the following table.

Step	Temperature and time	Cycles
1	50°C x 10 min	1
2	95°C x 3 min	1
3	95°C x 15 sec	40

	60°C x 10 sec	
	72°C x 15 sec	

Fluorescence should be collected at 60°C via Green, Yellow, Orange and Red channels.

11. Programming Other Machines

If you apply this kit to other Real-Time PCR machines, program it according to the following table (please note that following program is not applicable to RotorGene):

Step	Temperature and time	Cycles
1	50°C x 10 min	1
2	95°C x 3 min	1
3	95°C x 15 sec	40
	60°C x 30 sec	
	72°C x 15 sec	

Fluorescence should be collected at 60°C for FAM, VIC, ROX and Cy5 dyes.

12. Data Analysis: Rotor-Gene

To analyze data briefly, click on analysis menu and then under Quantitation tab double click on cycling A. Green. Manually put threshold at 0.1. Repeat the above for Yellow, Red and put threshold at 0.05 for Orange Channels.

Figures 1, 2, 3 and 4 represent typical graphs for Rotor-Gene.

To interpret the results, please note that:

- Signal in **Green** channel is due to **Influenza A**, **Yellow** channel due to **Influenza B**, **Red** channel due to **RSV** and **Orange** channel due to **IC or RNase P**.

Note that a sample is considered Positive only if it has a sigmoid graph and log phase, and only then CT is reliable and can be used.

In the absence of sigmoid graph and log phase, sample is considered Negative and CT if present is not reliable.

Consider following points when analyzing:

- A sample is **Positive** for Influenza A if it is positive in **Green** channel with sigmoid graphs and CT of less than 35 and also positive in **Orange** channel with a CT of 20-35.
- A sample is **Positive** for Influenza B if it is positive in **Yellow** channel with sigmoid graphs and CT of less than 35 and also positive in **Orange** channel with a CT of 20-35.
- A sample is **Positive** for RSV if it is positive in **Red** channel with sigmoid graphs and CT of less than 35 and also positive in **Orange** channel with a CT of 20-35.
- A sample is **Negative** for Influenza A, B and RSV if it is negative in Green and Yellow and Red channels while it is positive in Orange channel with a sigmoid graph and CT of 20-35.
- Results are **Inconclusive** and the test should be repeated if a sample is negative in all four channels.

Interpretation of results is summarized in table 1.

Green	Yellow	Red	Orange	Result
+	-	-	+	Positive for Flu A
-	+	-	+	Positive for Flu B
-	-	+	+	Positive for RSV
-	-	-	+	Negative
-	-	-	-	Inconclusive

Table 1. Interpretation of results


Fig 1. Typical Controls graph in Green channel for Rotor-Gene


Fig 2. Typical Controls graph in Yellow channel for Rotor-Gene


Fig 3. Typical Controls graph in Red channel for Rotor-Gene


Fig 4. Typical Controls graph in Orange channel for Rotor-Gene

13. Sensitivity

The analytical detection limit of the kit was assessed with cloned target of Influenza A, B and RSV genome and showed a limit of detection equal to 50 copy/ μ l for Influenza A, 250 copy/ μ l for Influenza B and 150 copy/ μ l for RSV.