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HCV Genotype RG Kit Manual

**For Genotyping of Hepatitis C Virus
(HCV) with Real-Time RT-PCR**

For use with Rotor-Gene
Research use only

NG-WI-ASL-03-102
Version 1.2
September 2017


Table of Contents:

1. Introduction	2
2. Kit Contents	3
3. Storage and Stability	4
4. Additionally Required Materials.....	4
5. General Precautions.....	4
6. Specimen, Storage and Transport.....	5
7. Interfering Substances	6
8. RNA Isolation	6
9. RT-PCR Protocol	6
10. Programming of Rotor-Gene.....	7
11. Programing Other Machines	7
12. Data Analysis	7
13. Sensitivity.....	10

HCV Genotype RG kit is intended for the genotyping of HCV type 1, 2 and 3 using RNA extracted from plasma. This kit is designed for use with Rotor-Gene and for research use only.

1. Introduction

Hepatitis C virus (HCV) is the major cause of cirrhosis and hepatocellular carcinoma. So far 6 major genotypes and more than 80 subtypes have been identified.

Based on current clinical trials, pre-treatment viral load and HCV genotype are two best indicators for medication choice as well as patient response to therapy. Therefore, once patient tested positive for HCV, both viral load and genotype should be determined. Moreover since HCV genotypes have significant differences in their global distribution and prevalence, genotyping can be a useful tool for epidemiologic studies.

The gold standard method for HCV genotyping is sequencing and phylogenetic studies. Since this method is not feasible for routine practice in a clinical laboratory, other methods have been used for this purpose which most of them are PCR-based.

HCV Genotype RG kit provides a ready-to-use One-Step Real-Time RT-PCR system for genotyping of HCV types 1, 2 and 3. In this method application of fluorescent probes allows detection of amplified product. Analysis of fluorescent kinetics also leads to quantitation and genotyping of the target sequence in the reaction without requiring post-amplification analysis, reducing the possibility of contamination with the PCR product. This kit is

designed for use with Rotor-Gene. HCV Genotype RG does not detect other HCV types including types 4, 5 and 6.

2. Kit Contents

The kit contains a manual, a CD with Rotor-Gene template and following reagents:

Label	Content	Quantity
HCV Gen Mix RG	RT-PCR Mix*	300 µl
HCV-1	HCV Type 1 control	150 µl
HCV-2	HCV Type 2 control	150 µl
HCV-3	HCV Type 3 control	150 µl
Water	PCR Grade Water	200 µl

* 1, 2 and 4 tubes for 24, 48 and 96 reaction kits.

3. Storage and Stability

The kit components should be shipped and stored at -20°C and are stable until the expiry date mentioned. Avoid repeated freeze-thaws more than few times to prevent reduced sensitivity.

4. Additionally Required Materials

To use this kit, you need the following items:

- Real-Time PCR machine and accessory computer
- Table top microtube centrifuge
- Vortex Mixer

- Dry Block Heater
- Adjustable pipettors and nuclease free filtered tips
- RNA extraction kit
- Nuclease free 1.7ml microtubes and PCR microtubes
- Disposable powder-free gloves
- Cold block

Note! This kit does not require cDNA synthesis reagents. They are already included in the Mix.

5. General Precautions

In order to prevent false results, always pay attention to the following points:

- **Treat all samples as potentially infectious.**
- Within the pre-PCR work area assign three separate spaces for: a) Sample storage and extraction; b) Reagent preparation where the HCV Gen Mix is aliquoted into tubes; and c) Reaction preparation area for addition of extracted RNA to the tubes.
- Always wipe the working surfaces with 70% Ethanol before and after work.
- **Thaw kit components on crushed ice** completely, mix by flickering followed by a quick spin and **store on crushed ice after.**
- Do not place PCR tubes on crushed ice. Use cold blocks instead.

6. Specimen, Storage and Transport

Peripheral blood should be collected in sterile condition in proper and sterile tubes. We recommend EDTA or citrate plasma for HCV detection and genotyping. Whole blood or plasma should be shipped at +4°C. Upon receipt, plasma should be separated from whole blood, then stored at +4°C for few days or aliquoted and stored at -20°C for up to few weeks.

7. Interfering Substances

Heparin (more than 10 IU/ml) affects the PCR. Blood collected in heparin containing tubes should not be used. Samples of heparinized patients must not be used as well. Elevated levels of bilirubin (≤ 4.5 mg/dl) and lipids (≤ 1000 mg/dl) and hemolytic samples do not influence the extraction and PCR.

8. RNA Isolation

RNA isolation can be performed using different kits from various manufacturers. We recommend using either of the following:

- High Pure Viral Nucleic Acid Kit (Cat. no. 11858874001, Roche Applied Science, Mannheim, Germany)
- QIAamp UltraSens ® Virus Kit (Cat. no. 53704, Qiagen GmbH, Hilden, Germany)
- QIAamp MiniElute Virus Spin Kit (Cat. no. 57704, Qiagen GmbH, Hilden, Germany)

9. RT-PCR Protocol

Thaw kit reagents on crushed ice completely followed by a brief mixing and a quick spin. **Pipette 12.5ul of HCV Gen Mix RG directly to each tube followed by adding 12.5ul of isolated RNA or Positive Ctrl or water.**

Cap the tubes and visually inspect to make sure all are capped securely. Place tubes in the machine.

Note: If using Rotor-Gene attach the locking ring too.

10. Programming Rotor-Gene

Before you start the machine make sure you have attached the locking ring on the rotor!

Open the CD provided in the kit and double click on “HCV Gen 0.2” or “HCV Gen 0.1” according to the used tubes. Program starts. Click on Start button (Green button on the top menu). On the pop up window click start again and save the run file.

11. Programming Other Machines

HCV Gen RG kit is intended for use with RotorGene machine. If you apply this kit to other Real-Time PCR machines, program it according to the following table:

Step	Temperature and time	Cycles
1	50°C x 30 min	1
2	95°C x 15 min	1
3	95°C x 15 sec	10
	62°C x 30 sec	
4	95°C x 15 sec	45
	60°C x 40 sec*	
	72°C x 15 sec	

* 15 seconds for Rotor-Gene machine.

Fluorescence should be collected at 60°C for FAM (Green), Hex (Yellow) and Cy5 (Red) during the final 45 cycles (section 4 of the profile).

12. Data Analysis

Analysis of data is relatively simple and quick. Each of HCV genotypes generates signal in a different channel (**Green channel: Type 1, Yellow channel: Type 3, Red channel: Type 2**).

Briefly, click on analysis menu and then under Quantitation tab double click on "Cycling A. Green". Close the pop up window for

Automatic Threshold and set threshold manually on 0.1 for Green channel and 0.2 for Yellow and Red channels.

Note that a sample is considered Positive only if it has a sigmoid graph and log phase, and only then CT is reliable and can be used.

In the absence of sigmoid graph and log phase, sample is considered Negative and CT if present is not reliable.

Consider following points when analyzing:

- A sample is **Positive for HCV Genotype 1** when it is positive in **Green channel**.
- A sample is **Positive for HCV Genotype 2** when it is positive in **Red channel**.
- A sample is **Positive for HCV Genotype 3** when it is positive in **Yellow channel**.

If a sample is Negative in all three Green, Yellow and Red channels, there are following possibilities:

- Sample is Negative for HCV.
- Sample is Positive for HCV genotypes 1, 2 or 3 but the viral load is lower than the kit sensitivity.
- Sample is Positive for HCV genotypes 4, 5 or 6.

Note! according to published studies, type 1 is the most prevalent type in Iran with type 3 as the second most prevalent type. Type 2 has also been reported. Types 4, 5 and 6 have not been reported from Iran yet.

Figures 1, 2 and 3 represent typical graphs for Rotor-Gene machine.


Figure 1: HCV Controls in Green Channel for Rotor-Gene


Figure 2: HCV Controls in Yellow Channel for Rotor-Gene


Figure 3: HCV Controls in Red Channel for Rotor-Gene

13. Sensitivity

The analytical detection limit of the kit was assessed with dilution series of the cloned HCV genotype 1 sequence and showed a limit of detection equal to 25 copy/ μ l.