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BRAF RQ

Kit Manual

**For Real-Time PCR
Detection of BRAF Mutations**

Research use only

NG-WI-ASL-40-100
Version 1.0
Fall 2020


Table of Contents:

1. Introduction 2

2. Kit Contents 2

3. Storage and Stability 3

4. General Precautions..... 3

5. Additionally Required Items 4

6. Specimen 4

7. DNA Extraction..... 4

8. Protocol: DNA Sample Assessment 5

9. Programming Real-Time PCR 5

10. Analysis: Sample Assessment..... 6

11. Protocol: Detection of BRAF mutation..... 8

12. BRAF Mutation Detection Analysis 9

13. Analytical Sensitivity16

BRAF RQ kit is intended for the detection BRAF V600 mutations in human genomic DNA. This kit is designed for use with Real-Time PCR machine. This kit is for research use only!

1- Introduction

BRAF oncogene is among the most frequently mutated kinases in human cancer. Mutations in codon V600 have been reported in different types of cancers including 40-50% of melanomas, 10-70% of thyroid carcinomas, 10% of colorectal cancers and 3-5% of non-small cell lung cancers (NSCLC).

Most of the BRAF mutations are located in codon 600 and constitute V600E, V600Ec, V600D, V600K, V600R.

BRAF RQ Kit provides ready to use reagents for detection of these five mutations based on Real-Time PCR technology. In this method application of fluorescent probes and analysis of fluorescent kinetics allows detection of amplified product without requiring post-amplification analysis. Therefore, test is accomplished much faster and possibility of contamination with the PCR product is reduced sharply. This kit also includes internal control for preventing false negative results due to failure in DNA extraction or PCR reaction.

2- Kit Contents

The kit contains a manual and CD and following reagents:

Label	Content	Quantity	Volume
BRAF Ctrl Mix	PCR Mix for quality control	2	480 µl
V600E Mix	PCR Mix to check V600E mutation	1	480 µl
V600EC Mix	PCR Mix to check V600Ec mutation	1	480 µl

V600D Mix	PCR Mix to check V600D mutation	1	480 μ l
V600K Mix	PCR Mix to check V600K mutation	1	480 μ l
V600R Mix	PCR Mix to check V600R mutation	1	480 μ l
BRAF Pos	1% Positive control	1	250 μ l
BRAF Neg	Negative control	1	250 μ l
Water	PCR Grade Water	1	200 μ l

3- Storage and Stability

The kit components should be shipped and stored at -20°C and are stable until the expiry date mentioned. Avoid repeated freeze-thaws especially for PCR Mix more than few times to prevent reduced sensitivity.

4- General Precautions

In order to prevent false results, always pay attention to the following points:

- **Treat all samples as potentially infectious.**
- Within the pre-PCR work area, assign three separate spaces for: a) Sample storage and extraction; b) Reagent preparation where the Master Mix is aliquoted into tubes; and c) Reaction preparation area for addition of templates to the tubes.
- Always wipe the working surfaces with 70% Ethanol before and after work.
- Thaw kit components on ice completely, mix by flickering followed by a quick spin and store on crushed ice after.
- Do not place PCR tubes on crushed ice. Use cold blocks instead.

5- Additionally Required Items

To use this kit, you need the following items:

- Real-Time PCR machine and accessory computer
- Adjustable pipettors and nuclease free filtered tips
- Vortex Mixer
- DNA extraction kit and required equipment/items
- PCR microtubes
- Disposable powder-free gloves
- Cold block

6- Specimen

Sections of fixed tissue (paraffin block) or unfixed tissue can be examined. Note that tumor tissue is heterogeneous with uneven distribution of mutant cells; therefore different parts of tissue may produce different results.

7- DNA Extraction

DNA extraction can be performed with different kits from various manufacturers. We recommend using:

- QIAamp Fast DNA Tissue Kit (Cat. no. 51404, Qiagen GmbH, Hilden, Germany)
- Roche; High Pure FFPE DNA Isolation Kit (Cat. No. 06 650 767 001, Roche Diagnostics GmbH, Mannheim, Germany)

Note that, the minimum volume of DNA required for the test is 40 microliters. We recommend extracting 100 microliters or higher volume of DNA, enough to repeat the test if necessary.

8- Protocol: DNA Sample Assessment

Before examining a sample for BRAF mutations, quality of DNA should be assessed. If the results are within the desired range, then the second test for detecting BRAF mutations will be performed. To qualify DNA extraction, follow steps below.

Note! Check the volume of extracted DNA sample and make sure it is more than 40 microliters and preferably about 100 microliters. Volumes less than 40 microliters are not enough to proceed.

First, thaw the **BRAF Ctrl Mix** on ice completely, mix by inversion followed by a quick spin and store on crushed ice after. Place the required number of microtubes on cold block including one for each sample, plus two for negative control and water.

Pipette 20µl of BRAF Ctrl Mix to each microtube. Continue by adding 5µl of DNA sample, Negative control and Water to each tube.

Cap the tubes and inspect visually to make sure all are capped securely. Place tubes in the machine.

9- Programming Real-time PCR

Rotor-Gene or StepOne could be setup using templates provided in the CD. Templates may also be downloaded from our website www.novingene.com. For Rotor-Gene, select “BRAF 0.2” or “BRAF 0.1” according to the PCR tubes.

You may also set up the machine as below.

Step	Temperature and time	Cycles
1	95°C x 3 min	1
2	95°C x 15 sec	45
	60°C x 60 sec	

Fluorescence should be collected at 60°C for FAM/Green and VIC/Yellow dyes/channels.

All mixes contain ROX with final concentration of 300nM in the reaction.

10-Analysis: Sample Assessment

Data Analysis is performed in two steps. First, the test validity is assessed and then DNA sample analysis. As the initial step, results of negative control will be evaluated.

Perform quantitative analysis for both **DNA sample (FAM/Green channel)** and **internal control (VIC/Yellow channel)** with the threshold set at 0.1 for both channels.

Note that, a sample is considered Positive only if it has a sigmoid graph and log phase, and only then CT is reliable and can be used. In the absence of sigmoid graph and log phase, sample can not be considered as positive (even CT present) result is not reliable.

Test is valid only if:

- 1) Water sample is **Negative** in FAM/Green channel and it is **Positive** in VIC/Yellow channel with CT of 27-33.
- 2) Negative controls is **Positive** in VIC/Yellow channel with CT of 27-33.
- 3) Negative control is **Positive** in FAM/Green channel with CT of 23-28.

Expected results for a valid test are summarized in Table1.

Sample	VIC/Yellow	FAM/Green
NTC	Pos (CT: 27-33)	Neg
Negative Ctrl	Pos (CT: 27-33)	Pos (CT: 23-28)

Table1. Expected results for valid Sample Assessment.

If above are met, then the test results are valid and may proceed to DNA sample assessment.

B) DNA sample Quality Analysis:

DNA sample is qualified only if sample is positive in FAM/Green channel with CT of 22-30 and positive in VIC/Yellow channel with CT of 27-33. This sample can be further examined for BRAF mutations.

Sample should be diluted with water, if a sample is positive in FAM/Green channel with CT of less than 22 and positive in VIC/Yellow channel with CT of 27-33. Dilute the sample to reach CT of 22-27. By 2X dilution of sample, CT is increased one cycle. After diluting the sample, proceed to the BRAF mutation test.

Result is invalid, if a sample in FAM/Green channel has CT above 30 or if VIC CT is below 27 or above 33. In this case, repeating DNA extraction is recommended.

DNA sample Quality Analysis is summarized in Table 2.

	FAM/Green	VIC/Yellow	Result
1	CT: 22-30	CT: 27-33	Valid
2	CT<22	CT: 27-33	Sample dilution
3	CT>30	CT<27 or CT>33	Invalid

Table2. DNA Quality Analysis

11-Protocol: Detection of BRAF mutations

To detect BRAF mutations, each sample should be examined with six mixes, one for DNA quantitation and one for each of five mutations. Therefore, to examine only one sample, 24 microtubes are required; six for the sample, 6 for the Positive control, six for the Negative control and six for NTC. Respectively for each extra sample, 6 microtubes will be added. Therefore, for two samples, 30 microtubes and for three samples 36 microtubes are required. Figure1 shows tube setup for three samples.

To start, Place required number of tubes on cold block organized in series of six each.

Pipette 20µl of BRAF Ctrl Mix to each tube in the first series.

Pipette 20µl of V600E Mix to each tube in the second series.

Pipette 20µl of V600Ec Mix to each tube in the third series.

Pipette 20µl of V600D Mix to each tube in the fourth series.

Pipette 20µl of V600K Mix to each tube in the fifth series.

Pipette 20µl of V600R Mix to each tube in the sixth series.

Continue by adding 5µl of extracted DNA, Positive control, Negative Control and water to each tube. Consider the first and second tube in each series for positive control and Negative control and the third for water/NTC. The next tubes would be for samples.

Cap the tubes and visually inspect to make sure all are capped securely. Place tubes in the machine and program it according to the section 9.

Note: if using Rotor-Gene attach the locking ring too.


Fig1. Tubes setup for Mixes and samples.

12-BRAF Mutation Detection Analysis

Before analyzing results for BRAF mutations, test validity should be verified. To do so, results of positive and negative controls and samples with Control Mix will be examined. If the test is valid, then may proceed to BRAF mutation analysis.

To analyze results, set threshold at 0.1 for FAM/Green and VIC/Yellow channel.

Note that a sample is considered Positive only if it has a sigmoid graph and log phase, and only then CT is reliable. In the absence of sigmoid graph and log phase, sample is considered Negative and CT if present is not reliable.

A) Test is valid only if:

- 1) NTC is Negative in FAM/Green channel with all Mixes, and
- 2) Positive Ctrl, Negative Ctrl and NTC are Positive in VIC/Yellow channel with each Mix with a CT of 27-33, and
- 3) Positive and Negative controls are Positive in FAM/Green channel with Control Mix with a CT of 22-27, and
- 4) Positive control is Positive in FAM/Green channel **with each of specific Mixes** with a CT of 30-35, and
- 5) Sample is Positive in VIC/Yellow channel with all mixes with CT of 27-33.
- 6) Sample is Positive in FAM/Green channel with Control Mix with a CT of 22-30 and

Above steps for Test validation are summarized in Figure 1.

If all above are met, then the Test is valid and results can be analyzed further for BRAF mutations.


Fig 1. Test validation flowchart

B) Data analysis for BRAF mutations

- 1) Select the patient's samples, which are positive in FAM/Green channel with BRAF specific mixes with a CT of 20-40, and document the CTs in the Table 3.
- 2) Calculate ΔCT for the above selected samples through the following equation.

$$\Delta CT = \text{Mutation Mix CT} - \text{Control Mix CT}$$

- 3) Document the ΔCT in the Table 3 too.
- 4) If the ΔCT is within the valid range as mentioned in the Table 3, the sample has the BRAF mutation.
- 5) If a sample is positive for two or more BRAF mutations, it is considered positive for the mutation with the lowest ΔCT and negative for the other mutations.

Above steps for detection of BRAF mutations are summarized in Figure 2.

BRAF Mix	Valid CT	Sample CT	Valid ΔCT	Sample ΔCT	Result
Ctrl Mix	22-30		-		Pass
600E Mix	20-40		≤ 7.3		V600E or V600K Pos
600Ec Mix	20-40		≤ 10.6		V600Ec Pos
600D Mix	20-40		≤ 9.6		V600D Pos
V600K Mix	20-40		≤ 9		V600K Pos
V600R Mix	20-40		≤ 8.3		V600R Pos

Table 3. Valid CT and ΔCT ranges for BRAF Mixes.

BRAF RQ (V1.0)

As an example, If CT of a sample is 24.7 with Control Mix, 32.2 with V600E Mix, 25.5 with V600Ec Mix, 37.4 with V600D Mix, 39.1 with V600K Mix and 43.5 for V600R mix, then ΔCT is 10.5 for V600E (35.2 - 24.7), 5.8 for V600Ec (30.5-24.7), 15.7 for V600D and 17.4 for V600K, 18.8 for V600R (Table. 4). According to table 3, delta CT for V600Ec is in valid range, and patient is Positive for V600Ec.

BRAF Mix	Valid CT	Sample CT	Valid ΔCT	Sample ΔCT	Result
Ctrl Mix	22-30	24.7	-		Valid
V600E Mix	20-40	35.2	≤ 7.3	10.5	Neg
V600Ec Mix	20-40	30.5	≤ 10.6	5.8	Positive
V600D Mix	20-40	40.4	≤ 9.6	15.7	Neg
V600K Mix	20-40	42.1	≤ 9	17.4	Neg
V600R Mix	20-40	43.5	≤ 8.3	18.8	Neg

Table 4. BRAF mutation analysis for a specific sample.

Note that a normal sample may cross-react with some of the mutation specific mixes. However, the calculated ΔCT will always fall outside the valid range. Therefore, a sample with CT of 20-40 in FAM/Green channel for BRAF mutation mixes, will only considered positive only if ΔCT is within the valid range.

Note: if sample is positive with V600E Mix, it could be positive for either of V600E and V600K Mutation.


Fig 2. Sample analysis flowchart for BRAF mutation detection

Results can also be interpreted as below

- Sample is **Negative** for BRAF mutations, if it is negative in FAM/Green channel with all mixes or has CT above 40.
- Sample is **Negative** for BRAF mutations if it is positive in FAM/Green channel with CT of 20-40 for one or more BRAF mixes, but ΔCT is higher than cutoff range (Table 3).

BRAF RQ (V1.0)

- Sample is **Positive** for BRAF mutations, if sample has a valid CT and ΔCT is lower than cutoff range (Table 3) in FAM/Green channel for one of the five BRAF mixes. Briefly,
 - It is Positive for V600E or V600K if with V600E Mix has a $CT \leq 40$ and $\Delta CT \leq 7.3$.
 - It is Positive for V600Ec if with V600Ec Mix has a $CT \leq 40$ and $\Delta CT \leq 10.6$.
 - It is Positive for V600D if with V600D Mix has a $CT \leq 40$ and $\Delta CT \leq 9.6$.
 - It is Positive for V600K if with V600K Mix has a $CT \leq 40$ and $\Delta CT \leq 9$.
 - It is Positive for V600R if with V600R Mix has a $CT \leq 40$ and $\Delta CT \leq 8.3$.
- If a sample has a valid CT and ΔCT (Table 3) for more than one Mix, sample is Positive for the BRAF mutation, which has the lowest ΔCT .

Interpretation of results for BRAF mutation test are summarized in Table 5.

Note that, if a sample is Negative with this kit in terms of BRAF mutations, the following conditions should be considered:

- Sample is Negative only for BRAF mutations which was mentioned and it could be Positive for other BRAF mutations.
- Sample is Positive for BRAF mutations in quantities below the sensitivity of the kit.

BRAF Mix	Sample CT	Sample Δ CT	Conclusion
600E Mix	>40	-	Neg for 600E mutation
	20-40	>7.3	Neg for 600E mutation
	20-40	$\leq$ 7.3	Pos for 600E or 600K mutation
600Ec Mix	>40	-	Neg for 600Ec mutation
	20-40	>10.6	Neg for 600Ec mutation
	20-40	$\leq$ 10.6	Pos for 600Ec mutation
600D Mix	>40	-	Neg for 600D mutation
	20-40	>9.6	Neg for 600D mutation
	20-40	$\leq$ 9.6	Pos for 600D mutation
V600K Mix	>40	-	Neg for 600K mutation
	20-40	>9	Neg for 600K mutation
	20-40	$\leq$ 9	Pos for 600K mutation
600R Mix	>40	-	Neg for 600R mutation
	20-40	>8.3	Neg for 600R mutation
	20-40	$\leq$ 8.3	Pos for 600R mutation

Table 5. Interpretation of results for BRAF mutation.

13-Analytical Sensitivity

Analytical sensitivity of this assay is equivalent to the percentage of mutated BRAF DNA that can be identified in the background of wild-type DNA and is mentioned in Table 6. The sensitivity depends on the DNA quantity in a sample. Maximum sensitivity of is achieved when a sample with Control Mix has a CT of 22-27 in FAM/Green channel; with CT of 28-30, the sensitivity will decrease for some mutations.

Reaction	Ctrl Mix CT: 22-27	Ctrl Mix CT: 28-30
V600E	1%	2%
V600Ec	0.5%	1%
V600D	0.5%	1%
V600K	1%	2%
V600R	2%	4%

Table 6. Sensitivity of the BRAF mutation detection assay.