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# **JAK2 RQ Kit Manual**


**For Real-Time PCR Detection of  
JAK2 V617F Mutation**

For use with Rotor-Gene or StepOne  
Research use only

NG-WI-ASL-10-301

Version 3.1

April 2018

**Table of Contents:**

|  |    |
|--|----|
| 1. Introduction .....                    | 2  |
| 2. Kit Contents .....                    | 3  |
| 3. Storage and Stability .....           | 3  |
| 4. General Precautions .....             | 3  |
| 5. Additionally Required Materials.....  | 4  |
| 6. Specimen, storage and transport ..... | 4  |
| 7. Interfering Substances .....          | 4  |
| 8. DNA isolation .....                   | 5  |
| 9. Internal control (IC).....            | 5  |
| 10. PCR Protocol .....                   | 5  |
| 11. Programming Rotor-Gene .....         | 6  |
| 12. Programming StepOne .....            | 6  |
| 13. Programming other machines .....     | 6  |
| 14. Data Analysis: Rotor-Gene .....      | 7  |
| 15. Data Analysis: StepOne .....         | 8  |
| 16. Sensitivity .....                    | 10 |

**JAK2 RQ** kit is intended for detection of JAK2 V617F mutation in human genomic DNA with Rotor-Gene 6000/Q or StepOne instrument. This kit is for research use only

## 1. Introduction

JAK2 (Janus Kinase 2) is a Tyrosine Kinase located in cytoplasm with essential role in signaling pathways for cytokines and growth factors. The acquired mutation G1849T replaces valine with phenylalanine (V617F). This substitution results in constitutively active JAK2 which leads to uncontrolled cell proliferation in the absence of growth factors. This mutation is found in the majority of BCR-ABL-negative myeloproliferative disorders (MPDs) and has become a main diagnostic test for polycythemia vera (PV), essential thrombocythemia (ET) and primary myelofibrosis (PMF). JAK2 RQ kit provides a ready-to-use Real-Time PCR system for detection of JAK2 V617F mutation with Rotor-Gene or StepOne machines. In this method application of fluorescent dye labeled probes allows detection of amplified product without requiring post-amplification analysis, reducing the possibility of contamination with the PCR product.

This kit also incorporates an *Internal Control* (IC) to evaluate DNA extraction quality or PCR inhibition.

## 2. Kit Contents

The kit contains a manual, a CD with Rotor-Gene and StepOne templates and following reagents:

| Label         | Content               | Quantity    |
|---------------|-----------------------|-------------|
| JAK2 Mix      | PCR Master mix*       | 480 $\mu$ l |
| Pos Ctrl 5%   | Positive Control 5%   | 50 $\mu$ l  |
| Pos Ctrl 0.1% | Positive Control 0.1% | 50 $\mu$ l  |
| Negative Ctrl | Negative Control      | 50 $\mu$ l  |
| Water         | PCR Grade Water       | 200 $\mu$ l |

\* 1, 2 and 4 tubes for 24, 48 and 96 reaction kits, respectively.

## 3. Storage and Stability

The kit components should be shipped and stored at -20°C and are stable until the expiry date mentioned. Avoid repeated freeze-thaws more than few times to prevent reduced sensitivity.

## 4. General Precautions

In order to prevent false results, always pay attention to the following points:

- **Treat all samples as potentially infectious.**
- Within the pre-PCR work area assign three separate spaces for: a) Sample storage and extraction, b) Reagent preparation where the Master Mix is aliquoted into tubes and c) Reaction preparation area for addition of extracted DNA to the tubes.
- Always wipe the working surfaces with 70% Ethanol before and after work.
- Thaw kit components on **crushed ice** completely, mix by flickering followed by a quick spin and **store on crushed ice while working.**

- Do not place PCR tubes on crushed ice. Use cooling blocks instead.

## 5. Additionally Required Materials

To use this kit, you need the following items:

- Real-Time PCR machine and accessory computer
- Table top microtube centrifuge
- Vortex Mixer
- Dry Block Heater
- Adjustable pipettors and nuclease free filtered tips
- DNA extraction kit
- Nuclease free 1.7ml microtubes and PCR microtubes
- Disposable powder-free gloves
- Cold block

## 6. Specimen, Storage and Transport

Peripheral blood is the suitable sample for JAK2 test. We recommend EDTA or citrate as anticoagulant.

DNA can directly be extracted from whole blood. Alternately to increase sensitivity buffy coat can be used. For optimum results, a sample should contain about 5 to 20 million WBC per milliliter.

Whole blood or buffy coat can be shipped and stored at +4°C for few days or can be aliquoted and stored at -20°C which is stable for few months.

## 7. Interfering Substances

Heparin (more than 10 IU/ml) affects the PCR. Blood collected in heparin containing tubes should not be used. Samples of heparinized patients must not be used as well.

Elevated levels of bilirubin ( $\leq 4.5$  mg/dl) and lipids ( $\leq 1000$  mg/dl) and hemolytic samples do not influence the extraction and PCR.

## 8. DNA Isolation

DNA isolation can be performed using different kits from various manufacturers. We recommend the following:

- High Pure PCR Template Preparation Kit (Cat# 11796828001, Roche Applied Science, Mannheim, Germany)
- QIAamp DNA Blood Mini Kit (Cat. no. 51104, Qiagen GmbH, Hilden, Germany)

To reach the sensitivity of 0.1%, extracted DNA should have a concentration of 25-50 ng/ul.

## 9. Internal Control (IC)

To examine DNA extraction quality as well as the presence of PCR inhibitors and to prevent false negative results, primers and probe for an *Internal Control* (a housekeeping gene) is included in the PCR Master Mix. Internal control should generate a CT of 22-25 in Rotor-Gene and a CT of 24-28 in StepOne. See analysis section for more details.

## 10. PCR Protocol

Thaw the reagents on ice completely followed by a brief mixing and a quick spin. Place required number of tubes on cooling or loading block. Consider one tube for each sample plus three for Positive control, Negative control and NTC.

**Pipette 20ul of JAK2 Mix directly to each tube followed by adding 5ul of control or isolated DNA.**

Cap the tubes and visually inspect to make sure all are capped securely. Place tubes in the machine.

*Note: Working with StepOne instrument, spin tubes briefly before loading on the block.*

*Note: If using Rotor-Gene attach the locking ring.*

## 11. Programming Rotor-Gene

*Before you start the machine, make sure you have attached the locking ring on the rotor!*

Open the CD provided in the kit and double click on "JAK2 0.2" or on "JAK2-strip" depending on the tubes used. Program starts. Click on Start button (Green button on the top menu). On the pop up window click start again and save program on desired location.

## 12. Programming StepOne

Open the StepOne software (V 2.\*). On the Set Up menu click on "Template" and select the file on CD provided with the kit. Click on "Plate Setup". Few samples along with controls are defined. You may change plate set up using right click options (copy, paste, clear). You may also add/remove samples or change sample names on "Define Targets and Samples" menu. When finished, click on "Start Run" and save the experiment on desired location. Instrument will start shortly.

## 13. Programming Other Machines

If you apply this kit to other Real-Time PCR machines, program it according to the following table:

| Step | Temperature and time | Cycles |
|------|----------------------|--------|
| 1    | <b>95°C x 10 min</b> | 1      |
| 2    | <b>95°C x 15 sec</b> | 50     |
|      | <b>60°C x 60 sec</b> |        |

Fluorescence should be collected at 60°C for FAM and VIC/HEX dyes.

JAK2 Mix contains 300nM ROX in final reaction.


### 14. Data Analysis: Rotor-Gene

Analyze data according to manufacturer recommendations. Perform qualitative analysis for **JAK2 (Green channel)** and **Internal Control (Yellow channel)**. Briefly, click on "Analysis" menu and then under "Quantitation" tab, double click on "Cycling A. Green".


Close the pop up for Automatic Threshold and set the threshold on 0.1 and above background or Negative control.

Repeat the above for "Cycling A. Yellow" and manually set threshold on 0.1.

Refer to figures 1 and 2 for typical graphs.

**Figure 1.** Typical graph of JAK2 controls in green channel for Rotor-Gene

**Figure 2.** Typical IC graph in yellow channel for Rotor-Gene

Consider following points when analyzing:

**Note that a sample is considered Positive only if it has a sigmoid graph and log phase, and only then CT is reliable and can be used.**


**In the absence of sigmoid graph and log phase, sample is considered Negative and CT if present is not reliable.**

- A sample is **Positive** when it is positive in Green channel with a sigmoid graph and also positive in Yellow channel with a sigmoid graph and CT of 22-25.
- A sample is **Negative** if it is negative in Green channel but, positive in Yellow channel with a sigmoid graph and CT of 22 to 25.
- Results are **inconclusive** if a sample is negative in both of Green and Yellow channels. Improper DNA extraction or test set up could be the cause.

- Results are **inconclusive** if a sample is negative in Green channel but, positive in Yellow channel with CTs higher than 25. Low DNA concentration could be the cause.


### 15. Data Analysis: StepOne

Analyze data according to manufacturer recommendations. Briefly, click on "Analyze" and set the threshold at 0.1 for **JAK2 (FAM channel)** and 0.05 for **Internal Control (VIC channel)**. Refer to figures 3 and 4 for typical graphs.

**Figure 3.** Typical graph for JAK2 controls in FAM Channel for StepOne

**Figure 4.** Typical IC graph in VIC channel for StepOne

**Note that a sample is considered Positive only if it has a sigmoid graph and log phase, and only then CT is reliable and can be used.**

**In the absence of sigmoid graph and log phase, sample is considered Negative and CT if present is not reliable.**

Consider following points when analyzing:

- A sample is **Positive** when it is positive in FAM channel with a sigmoid graph and also positive in VIC channel with a sigmoid graph and CT of 24-28.
- A sample is **Negative** if it is negative in FAM channel but, positive in VIC channel with a sigmoid graph and CT of 24-28.
- Results are **inconclusive** if a sample is negative in both of FAM and VIC channels. Improper DNA extraction or test set up could be the cause.

- Results are **inconclusive** if a sample is negative in FAM channel but, positive in VIC channel with CTs higher than 28. Low DNA concentration could be the cause.

## 16. Sensitivity

The analytical detection limit of the kit was assessed with dilution series of the positive target in negative DNA and showed a limit of detection equal to 0.1% or 1 in 1,000. In other words, if in 1,000 cells only one is positive, it will be detected in 95% cases. To meet this level of sensitivity, extracted sample should contain 25-50 ng/ $\mu$ l DNA.