



Gibco™

Fetal Bovine Serum, Premium Plus



Catalog number: A4766801

Gibco Premium Plus Fetal Bovine Serum (FBS) is our highest quality FBS, meeting selective specifications to help provide consistent, reproducible results. Only Premium Plus FBS delivers our best quality specifications and testing levels.

Gibco sera: trusted quality for consistent performance

Premium sera features:

- Use with a broad range of cell types, especially sensitive cell lines
- Our most characterized sera with the lowest endotoxin and hemoglobin levels
- Endotoxin level: ≤ 1 EU/mL
- Hemoglobin level: ≤ 15 mg/dL
- Origin: United States, Australia, or New Zealand (refer to CoA for specific origin)

Additional features of this serum:

- Meets USP/EP guidelines with up to 96 harmonized quality specification tests
- Additional proliferation tests on top six cell lines
- EMA virus testing upon request

- Lowest risk of BSE and lower viral risk
- ISO 13485-certified; processed in FDA-registered facilities
- Triple-filtered at 0.1 μm
- Gibco bottle is easy to use in a hood, helps reduce risk of contamination, and allows you to perform cell culture more consistently
- Peel-off sticker on label provides handy reference for lot number and expiration date recording in lab notebook

Specifications	
Certifications/Compliance	Meets USP/EP guidelines; ISO 13485 Certified
Classification	Premium
Age	Fetal
Cell Type	Use for a broad range of cell types, especially sensitive cell lines
Country of Origin	Refer to CoA for specific origin
Form	Liquid
Format	Bottle
Grade	Premium Plus
Origin	United States, Australia Or New Zealand (Refer to CoA for Specific Origin)
Product Type	Fetal Bovine Serum
Quantity	500 mL
Size	500 mL
Species	Cattle/Bovine

Sterility	Sterile
Sterilization Method	Triple Filtered, 0.1 µm
Tested For	Meets USP/EP guidelines with up to 96 harmonized quality specification tests; Additional proliferation tests on top 6 cell lines; EMA virus testing upon request
Endotoxin Level	≤1 EU/mL
Green Features	Less Waste
Hemoglobin Concentration	≤15 mg/dL
Hemoglobin Level	≤15 mg/dL
Product Line	Gibco™
Purity or Quality Grade	Premium Plus
Shelf Life	5 years
Shipping Condition	Frozen
Storage Requirements	≤-10°C
Storage	≤-10°C

Contents & Storage

Store at or below -10°C.