

# **SAFETY DATA SHEET**

[Required under safety and health regulations for shipping and handling]

Version: 2019

Date Updated: June 24, 2019

#### SECTION 1. - - - - - PRODUCT AND COMPANY IDENTIFICATION - - - - - - -

**Product Name** 1,10-Phenanthroline, Monohydrate

Product Code(s) PB0693

Recommended Use For Laboratory Research Use Only

Not for Human or Animal Drug Use

**Supplier** Bio Basic Inc.

Address 20 Konrad Crescent, Markham, Ontario,

Canada, L3R 8T4

 Telephone
 (905) 474 4493

 Fax
 (905) 474 5794

 For Chemical Emergency Phone#
 (416) 995 9730

#### SECTION 2. ----- HAZARDS IDENTIFICATION -----

## **Emergency Overview**

#### **WHMIS Classification**

D1B Toxic Material Causing Immediate and Serious Toxic

Toxic by ingestion

**Toxic Effects** 

#### **GHS Classification**

Acute toxicity, Oral (Category 3)
Acute aquatic toxicity (Category 1)
Chronic aquatic toxicity (Category 1)

#### GHS Label elements, including precautionary statements

Pictogram

Signal word Danger

Hazard statement(s)

H301 Toxic if swallowed.

H410 Very toxic to aquatic life with long lasting effects.

Precautionary statement(s)

P273 Avoid release to the environment.

P301 + P310 IF SWALLOWED: Immediately call a POISON CENTER or doctor/ physician. P501 Dispose of contents/ container to an approved waste disposal plant.

**HMIS Classification** 

Health hazard: 2 Flammability: 0 Physical hazards: 0

#### **Potential Health Effects**

**Inhalation** May be harmful if inhaled. May cause respiratory tract irritation.

**Skin** May be harmful if absorbed through skin. May cause skin irritation.

**Eyes** May cause eye irritation. **Ingestion** Toxic if swallowed.

#### SECTION 3. - - - - COMPOSITION/INFORMATION ON INGREDIENTS - - - -

Chemical Name	EC No.	CAS-No	Weight %
1,10-Phenanthroline, Monohydrate	200-629-2	5144-89-8	<100

#### SECTION 4. - - - - - FIRST-AID MEASURES- - - - - -

#### General advice

Consult a physician. Show this safety data sheet to the doctor in attendance. Move out of dangerous area.

#### If inhaled

If breathed in, move person into fresh air. If not breathing, give artificial respiration. Consult a physician.

#### In case of skin contact

Wash off with soap and plenty of water. Take victim immediately to hospital. Consult a physician.

#### In case of eye contact

Flush eyes with water as a precaution.

#### If swallowed

Never give anything by mouth to an unconscious person. Rinse mouth with water. Consult a physician.

#### SECTION 5. ----- FIRE FIGHTING MEASURES -----

#### Conditions of flammability

Not flammable or combustible.

## Suitable extinguishing media

Use water spray, alcohol-resistant foam, dry chemical or carbon dioxide.

#### Special protective equipment for firefighters

Wear self-contained breathing apparatus for firefighting if necessary.

#### **Hazardous combustion products**

Hazardous decomposition products formed under fire conditions. - Carbon oxides, Nitrogen oxides (NOx)

#### Explosion data - sensitivity to mechanical impact

No data available

### Explosion data - sensitivity to static discharge

No data available

## SECTION 6. - - - - - ACCIDENTAL RELEASE MEASURES - - - - - -

#### Personal precautions

Wear respiratory protection. Avoid dust formation. Avoid breathing vapours, mist or gas. Ensure adequate ventilation. Evacuate personnel to safe areas. Avoid breathing dust.

## **Environmental precautions**

Prevent further leakage or spillage if safe to do so. Do not let product enter drains. Discharge into the environment must be avoided.

#### Methods and materials for containment and cleaning up

Pick up and arrange disposal without creating dust. Sweep up and shovel. Keep in suitable, closed containers for disposal.

#### SECTION 7. ----- HANDLING AND STORAGE-----

#### Precautions for safe handling

Avoid contact with skin and eyes. Avoid formation of dust and aerosols.

Provide appropriate exhaust ventilation at places where dust is formed.

#### Conditions for safe storage

Keep container tightly closed in a dry and well-ventilated place.

Hygroscopic. Store under inert gas. Keep in a dry place.

#### SECTION 8. - - - - EXPOSURE CONTROLS/PERSONAL PROTECTION- - - -

#### Personal protective equipment

#### Respiratory protection

Where risk assessment shows air-purifying respirators are appropriate use a full-face particle respirator type N99 (US) or type P2 (EN 143) respirator cartridges as a backup to engineering controls. If the respirator is the sole means of protection, use a full-face supplied air respirator. Use respirators and components tested and approved under appropriate government standards such as NIOSH (US) or CEN (EU).

## **Hand protection**

Handle with gloves. Gloves must be inspected prior to use. Use proper glove removal technique (without touching glove's outer surface) to avoid skin contact with this product. Dispose of contaminated gloves after use in accordance with applicable laws and good laboratory practices. Wash and dry hands.

#### Eye protection

Face shield and safety glasses Use equipment for eye protection tested and approved under appropriate government standards such as NIOSH (US) or EN 166(EU).

#### Skin and body protection

Complete suit protecting against chemicals, The type of protective equipment must be selected according to the concentration and amount of the dangerous substance at the specific workplace.

## Hygiene measures

Avoid contact with skin, eyes and clothing. Wash hands before breaks and immediately after handling the product.

#### Specific engineering controls

Use mechanical exhaust or laboratory fumehood to avoid exposure.

## SECTION 9. ----- PHYSICAL AND CHEMICAL PROPERTIES -----

## **Appearance**

Form powder Colour white

#### Safety data

pH No data available

Melting point/range: 100 - 104 °C (212 - 219 °F) - lit.

point/freezing point

Boiling point No data available
Flash point Not applicable
Ignition temperature No data available
Auto-ignition No data available

temperature

Lower explosion limit No data available
Upper explosion limit No data available
Vapour pressure No data available

Density No data available

Water solubility soluble

Partition coefficient: n-octanol/water

No data available

Dalada

Relative vapour

No data available

density

Odour No data available
Odour Threshold No data available
Evaporation rate No data available

#### SECTION 10. ------STABILITY AND REACTIVITY -----

#### **Chemical stability**

Stable under recommended storage conditions.

### Possibility of hazardous reactions

No data available

#### Conditions to avoid

Exposure to moisture

#### Materials to avoid

Acids, Strong oxidizing agents

#### **Hazardous decomposition products**

Hazardous decomposition products formed under fire conditions. - Carbon oxides, Nitrogen oxides (NOx) Other decomposition products - No data available

## SECTION 11. ----- TOXICOLOGICAL INFORMATION -----

## **Acute toxicity**

#### Oral LD50

LD50 Oral - Rat - 132 mg/kg

#### Inhalation LC50

No data available

#### **Dermal LD50**

No data available

#### Other information on acute toxicity

No data available

#### Skin corrosion/irritation

No data available

#### Serious eye damage/eye irritation

No data available

#### Respiratory or skin sensitisation

No data available

#### Germ cell mutagenicity

No data available

### Carcinogenicity

IARC: No component of this product present at levels greater than or equal to 0.1% is identified as

probable, possible or confirmed human carcinogen by IARC.

ACGIH: No component of this product present at levels greater than or equal to 0.1% is identified as a

carcinogen or potential carcinogen by ACGIH.

### Reproductive toxicity

No data available

## **Teratogenicity**

No data available

## Specific target organ toxicity - single exposure (Globally Harmonized System)

No data available

## Specific target organ toxicity - repeated exposure (Globally Harmonized System)

No data available

## **Aspiration hazard**

No data available

#### Potential health effects

**Inhalation** May be harmful if inhaled. May cause respiratory tract irritation.

**Ingestion** Toxic if swallowed.

**Skin** May be harmful if absorbed through skin. May cause skin irritation.

**Eyes** May cause eye irritation.

## Synergistic effects

No data available

## Additional Information

RTECS: Not available

#### SECTION 12. ----- ECOLOGICAL INFORMATION -----

#### **Toxicity**

No data available

### Persistence and degradability

No data available

## **Bioaccumulative potential**

No data available

#### Mobility in soil

No data available

#### PBT and vPvB assessment

No data available

#### Other adverse effects

An environmental hazard cannot be excluded in the event of unprofessional handling or disposal.

Very toxic to aquatic life with long lasting effects.

#### SECTION 13. ----- DISPOSAL CONSIDERATIONS -----

#### **Product**

Offer surplus and non-recyclable solutions to a licensed disposal company. Contact a licensed professional waste disposal service to dispose of this material. Dissolve or mix the material with a combustible solvent and burn in a chemical incinerator equipped with an afterburner and scrubber.

## Contaminated packaging

Dispose of as unused product.

## SECTION 14. ----- TRANSPORT INFORMATION -----

#### TDG (Canada)

UN number: 2811 Class: 6.1 Packing group: III

Proper shipping name: Toxic solids, organic, n.o.s. (o-Phenanthroline hydrate)

Reportable Quantity (RQ): Marine pollutant: No

Poison Inhalation Hazard: No

#### **IMDG**

UN number: 2811 Class: 6.1 Packing group: III EMS-No: F-A, S-A

Proper shipping name: TOXIC SOLID, ORGANIC, N.O.S. (o-Phenanthroline hydrate)

Marine pollutant: No

#### IATA

UN number: 2811 Class: 6.1 Packing group: III

Proper shipping name: Toxic solid, organic, n.o.s. (o-Phenanthroline hydrate)

## SECTION 15. ----- REGULATORY INFORMATION -----

#### **WHMIS Classification**

D1B Toxic Material Causing Immediate and Serious Toxic by ingestion

**Toxic Effects** 

This product has been classified in accordance with the hazard criteria of the Controlled Products Regulations and the MSDS contains all the information required by the Controlled Products Regulations.

#### SECTION 16. ----- OTHER INFORMATION-----

Further information: no limited for paper copy, just for internal uses.

For research use only. Not intended for human or animal diagnostic or therapeutic uses.

#### Disclaimer

The information provided on this SDS is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guide for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered as a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other material or in any process, unless specified in the text.

**End of SDS**