

AUTOCLAVABLE



## AUTOCLAVABLE ADJUSTABLE VALUM PIPETTE

This pipette is a general purpose pipette for sampling and dispensing accurate liquid volumes. The pipettes operate on the air displacement principle and use disposable tips. The pipette is autoclavable

( 121°C , 0.15 Mpa )

The pipettes cover a volume range of 0.5 µl to 10 ml .

All pipettes have been quality tested according to our company's demands. The quality control according to our company's involves gravimetric testing of each pipette with distilled water at 22°C using the manufacturer's original tips.

Cat. No.	Volume Range (ul)	Increment (ul)	Test Volume (ul)	Inacc (%)	Impr(%)	Tip
FA3013	0.5-10µl	0.1	10µl	±3.0	≤2.0	10µl
			5µl	±6.0	≤3.0	
			1µl	±10.0	≤3.0	
FA3014	2-20µl	0.5	20µl	±3.0	≤1.5	200µl
			10µl	±3.0	≤2.0	
			2µl	±10	≤3.0	
FA3015	5-50µl	0.5	50µl	±2.4	≤1.5	200µl
			25µl	±2.8	≤2.0	
			5µl	±6.0	≤3.0	
FA3016	10-100µl	1	100µl	±1.5	≤1.0	200µl
			50µl	±2.2	≤1.5	
			10µl	±3.0	≤2.0	
FA3017	20-200µl	1	200µl	±1.5	≤1.0	200µl
			100µl	±1.5	≤1.0	
			20µl	±3.0	≤1.5	
FA3018	50-200µl	1	200µl	±1.5	≤1.0	200µl
			100µl	±1.5	≤1.0	
			50µl	±2.2	≤1.5	
FA3019	100-1000µl	5	1000µl	±1.2	≤0.5	1000µl
			500µl	±1.5	≤0.5	
			100µl	±1.5	≤1.0	
FA3020	200-1000µl	5	1000µl	±1.2	≤0.5	1000µl
			500µl	±1.5	≤0.5	
			200µl	±1.5	≤1.0	
FA3021	1000-5000µl	50	5000µl	±1.0	≤0.5	5000µl
			2500µl	±1.0	≤0.5	
			1000µl	±1.2	≤0.5	
FA3022	2000-10000µl	0.1ml	10000µl	±0.5	≤0.2	10ml
			5000µl	±0.8	≤0.3	
			2000µl	±1.0	≤0.3	

Cat. No	Volume (ul)	Inacc. %	Impr.%
FA3001	5	±3.8	≤3.0
FA3002	10	±3.0	≤2.0
FA3003	20	±2.4	≤1.0
FA3004	25	±2.4	≤1.0
FA3005	50	±2.4	≤1.0
FA3006	100	±1.5	≤0.8
FA3007	200	±1.2	≤0.7
FA3008	250	±1.2	≤0.7
FA3009	500	±1.2	≤0.7
FA3010	1000	±1.0	≤0.5
FA3011	2000	±1.0	≤0.5
FA3012	5000	±1.0	≤0.5