

Multichannel pipette ME and MT



Multichannel pipette ME/MT

The fully autoclavable PIPETTES ME/MT are specially designed for working with gloves under sterile conditions. The volume range from 1µl to 300 µl is covered by 4 models each in the 8- and 12-channel versions.

The patented manifold system ensures individual tight sealing of the pipette tips on each channel.

- perfect measuring results
- fully autoclavable and UV-light resistant
- volume settings by pushbutton avoids snatching the gloves
- ergonomic design, handle can be turned to an individual convenient position
- patented system makes manual corrections of the tip sealing redundant
- maintenance free: all components are manufactured either from stainless steel or chemical resistant plastics
- each pipette comes with quality certificate, holder and calibration key

8-channel-pipette ME

Volume range µl	Model	Accuracy max. vol %	Precision < %	Quantity per Pack	Order No.
1-10	ME10	2,0	1,2	1	47 70 064
5-50	ME50	1,6	0,6	1	47 70 065
20-200	ME200	1,0	0,6	1	47 70 066
50-300	ME300	1,0	0,6	1	47 70 067

12-channel-pipette MT

Volume range µl	Model	Accuracy max. vol %	Precision < %	Quantity per Pack	Order No.
1-10	MT10	2,0	1,2	1	47 70 068
5-50	MT50	1,6	0,6	1	47 70 069
20-200	MT200	1,0	0,6	1	47 70 070
50-300	MT300	1,0	0,6	1	47 70 071